



MANGAR Pillowlifts

The Handy Pillowlift and the Sit-u-Up assist people who need help to rise from lying to sitting.

They minimise moving and handling risks and reduce the risk of musculoskeletal injury for the carer. With no fixtures or fittings, they sit on top of the mattress and can be used on any size bed, so it causes no disturbance to a second person in a double bed. Inflation and deflation can be stopped at any time and can lift a person on their back or side, allowing freedom to choose the most comfortable position.

TECHNICAL SPECIFICATIONS

Product	Handy Pillowlift	Sit -u-Up
Max user weight	135 kg	150 kg
Min height	Flat	Flat
Total length	80 cm	75 cm
Width	61 cm	70 cm
Total product weight	2 kg	4.8 kg
Back/headrest length	73 cm	73 cm
Backrest angle adjustment	Flat to 55°	Flat to 65°
Airflo Compressor	12	12
Lifting performance per charge	20 lifts	20 lifts



Features

1. Adjustable back support.
2. Easy to wash and change covers.
3. Spreader bar positioned against headboard to help secure Pillowlift in place.
4. Anti-slip pillow holder.
5. Push button hand control to inflate and deflate Pillowlifts.
6. Made from durable polyurethane coated nylon material.




Accessories

- **Pillowlift Quilt**
- **Pillowlift Cover**



Product codes:
Handy Pillowlift: MPCA120500
Sit-u-Up: MPCA130500



For more information visit
www.mangarhealth.com
 Email: enquiries@winnicare.uk
 Phone: 0800 2800 485

