



The Luna hoist is the only one that seamlessly fits onto any competitor track—saving you time, money, and disruption

## Luna, shaping the future of care

Upgrading, repairing and servicing your ceiling hoist shouldn't mean costly downtime or a full track replacement. Our innovative Luna hoist is the only one that seamlessly fits onto any competitor track—saving you time, money, and disruption. Keep patient care uninterrupted while benefiting from a quick and hassle-free upgrade.

Like any frequently used equipment, as hoists come to the end of their useable life, they will become increasingly vulnerable to costly breakdowns. Whether it's a repair or spare part that's required to get the hoist working again, there are often delays to a resolution, which will impact continuity of care for both patients and their cares.

For a quick, efficient solution, that minimises downtime and adds years of life to your CLS, use the Luna. They can be quickly and easily fitted onto any existing tracks, reducing delays in care that can be experienced waiting for spare parts and further engineer time.


We keep the Luna and the carts in stock in the UK, they can be shipped same day and delivered next day enabling your customers to keep their ceiling lift systems available, patients moving and care team efficient.



## A minimal design for a maximum lifting height

Best performance on lifting height and strap length

Fits in every environment, including with low ceiling height!



185 mm max.

Low height tracks

System height 241 mm

Patented Bayonett design

Innovative spreader bar design

Installing a Luna means immediate functionality and efficiency, leading to enhanced patient care and operational continuity.

- ✓ The Luna is beautifully designed, can be used with any in-situ track and is **simple to use**.
- ✓ The patented bayonet system and high reach strap (2500 mm) give an **unrivalled lifting height**.
- ✓ Carts are **extremely light and compact**, weighing only 5.2 kg.
- ✓ Enables a **fixed or portable lift**, delivering a 150 kg, 200 kg or 275 kg lifting load.
- ✓ **Fast to install and service** e.g. change of a battery, strap, hand control cord or clip, can all be done in less than 5 minutes and with 1 screwdriver
- ✓ Smart maintenance – replace only what needs to be changed – **60% less maintenance costs** than other systems.

For further information on the Winncare product range, please email [enquiries@winncare.uk](mailto:enquiries@winncare.uk)